

# Tissue&cell Genomic DNA Midi Purification Kit

## Description:

The Tissue & Cell Genomic DNA Midi Purification Kit provides a rapid, simple and effective approach to isolate the genomic DNA from various animal tissues up to 100mg and culture cells up to  $2 \times 10^7$ . The process is based on a spin column format; the procedure involves cell lysis with proteinase K digestion, nucleic acids absorption and DNA elution. There is no requirement for phenol/chloroform extraction or alcohol precipitation. Typical yield ranges from 100-120µg. DNA purified with this kit is suitable for various applications, including PCR, restriction enzyme digestion, cloning, dot blot analysis, etc.

## Characteristics:

- Fast spin-column
- Centrifugation-based method.
- No phenol/chloroform extraction and ethanol precipitation required.
- Consistent, high yields
- Complete removal of contaminants and inhibitors
- Operation time: < 60 minutes
- Wide range sample from plant tissue

Product Name	Size	Cat.No
Tissue&cell Genomic DNA Midi Purification Kit	25 Reactions	Bio-GTD-25



## Applications

- PCR
- Southern Blotting
- Medicolegal Analysis
- Real-Time PCR,
- AFLP, RFLP, PAPD.

## Storage Condition:

store at RT

## Components of the kit :

	<b>GTD-25</b>
1. Lysis Buffer	30 ml
2. Binding Buffer	30 ml
3. Proteinase K solution	10 mg/mL , 2.7mL
4. Wash Buffer I	30 ml (add 165ml of Ethanol before use)
5. Wash Buffer II	30 ml
6. Elution Buffer	30 ml
7. Midi Spin column	25 pcs
8. Midi Collection tube	25 pcs

### **Before beginning this procedure .:**

- Store the Proteinase K at -20°C . store all other kit components at RT
- RNA-free genomic DNA is required, prepared 100mg/ml RNase A ( **not include in the kit** ) with H<sub>2</sub>O.
- Turn on the water baths or heat blocks at 56°C and 70°C .

### **General Procedure :**

#### **Materials to be supplied by the user:**

- For tissue grinding : Small homogenizer ( fisher Tissue Tearor. Polytron or Turax). Alternatively, mortar and pestle.
- DNase-free RNaseA
- Trypsin (for adherent tissue culture cells only).
- PBS buffer or TE buffer ( culture cells only )
- Ethanol

## **Protocol for isolation of genomic DNA from tissue & Cell**

### **1. a. Tissue culture cells : ( < 2 X 10<sup>7</sup> )**

**\* Do not use too much sample, which will clog the column and obtain lower yield and quality.**

( I ) Harvest the cells (for adherent cells, trypsinize the cells before harvesting) and transfer to 50ml centrifuge tube.

( II ) Centrifuge at 15,000g for 10s to pellet the cells, remove the supernatant.

( III ) Add 1ml of PBS or TE buffer to the pellet, vortex or pipet to completely resuspend the pellet.

**Optional:** If RNA-free DNA is desired, add 20 µl RNaseA solution and incubate at RT for 5 min.

( IV ) Add 100 µl of **Proteinase K Solution** and 1ml **Binding Buffer**(vortex before use), vortex thoroughly to lyse the cells and proceed **directly** to **Step 6**.

### **b. Animal tissue ( <100 mg, < 80 mg for spleen) or insect ( < 100 mg ):**

**\*Do not use too much sample, which will clog the column and obtain lower yield and quality.**

Three methods can be used for tissue treatment:

( I ) Cut the tissue into small pieces, add 1 ml of **Lysis Buffer**, Proceed to **Step 2**.

( II ) Add 1 ml of **Lysis Buffer** to tissue and homogenize for 10s using a small homogenizer, then transfer to 50ml centrifuge tube. Proceed to **Step 2**.

( III ) Tissue may be ground in liquid nitrogen using a mortar and pestle. After grinding, transfer the ground tissue to 50ml centrifuge tube, then add 1 ml of **Lysis Buffer** and proceed to **Step 2**.

### **2. Add 100 µl of **proteinase K stock solution**, mix by vortexing.**

*\* Do not premix Lysis Solution and proteinase K solution before use, proteinase K may undergo self-digestion without substrate.*

3. Incubate at 56°C in water-bath or incubator for 1- 3 h or longer **until completely lysed.**, pulse vortex 5-10 sec occasionally during incubation.

4. **Optional:** If RNA-free DNA is desired, add 20 µl RNaseA solution and incubate at RT for 5 min.

5. Add 1 ml of **Binding Buffer**(vortex before use) into reaction, mix well by vortex.

6. **Incubating at 70°C in water-bath or heating block for 15 min.**

*\* A white precipitate may form after addition of Binding Buffer, which will redissolve during incubation at 70°C, and will not affect the DNA binding.*

*\* Some tissue debris ( i.e. zebra fish bone or insect exoskeleton ) may not be digestible, it is important to remove the debris by centrifuging at (10-15,000 ×g) for 5 min before loading on column, since these debris will clot the column.*

7. Add 1 ml **Ethanol**, mix thoroughly by vortex and apply the solution to midi spin column with collection tube, spin at top speed (10-15,000 ×g) for 5 min.

*\* A precipitate may form after addition of Ethanol, **apply all the solution and precipitate** to the column.*

8. Discard the flow-through, wash **twice** with 3 ml of **Wash Buffer** by spin at 15,000 ×g for 5 min.

9. Discard the flow-through, wash with 1ml of **Wash Buffer II** by centrifuge at 10,000×g for 5 min. to remove ethanol.

10. Discard the flow-through, Elution twice with 250-500uL **Elution Solution** or **H<sub>2</sub>O** ( pH must between 7.0-8.5 ), Incubate at 56°C for 5 min, to elute the DNA by centrifugation for 5 min, and store the DNA at -20°C.

*\* When using water to elute, make sure the pH value is within 7.0-8.5. Lower pH may cause lower DNA recovery.*