

PRODUCT INFORMATION

Product name : 100 mM dNTPs (Deoxynucleotide Triphosphates)
(high concentration)

Catalog number: **Bio-NTP-100 (A, C, G, T)**

Size : 4 separate vials containing 40 umoles of each dNTP in a volume of 400 ul

Description : Separate vials of dATP, dGTP, dCTP and dTTP each at a concentration of 100 mM. The dNTPs are made by dissolving ultra-pure triphosphate powders in buffer at pH 7.0

Storage : The dNTPs are stable for 6 months at -20°C, or long term storage at -70°C.

Application : PCR applications, DNA labelling and DNA sequencing.



BioKit Biotechnology Incorporation
Tel : 886-3-7668266
Fax : 886-3-7668192
E-mail: biokit@biokit.com.tw